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en INSTRUCTIONS FOR USE
Bipolar cables and adapter for electrosurgery

REF. No.	Generator type	Instrument type	Max. voltage [V _p]
Cables			
101-000	Erbe VIO / ICC / ACC / T series	EU forceps, EU flat connector	550
101-040	Erbe VIO / ICC / ACC / T series	BOWA forceps, EU flat connector	550
101-045	Erbe VIO / ICC / ACC / T series	US forceps	550
101-140	BOWA	BOWA forceps, EU flat connector	550
101-245	Erbe VIO *) / ICC / ACC / T series	ErgoLAP bipolar	600
101-445	Erbe VIO / ICC / ACC	Aesclapal forceps, flat connector	
101-345	Erbe VIO / ICC / ACC	Storz bipolar laparoscopic handle	550
287-040	Martin	BOWA forceps, EU flat connector	550
287-045	Martin	US forceps	550
287-050	Martin	BOWA forceps, EU flat connector	550
287-245	Martin	ErgoLAP bipolar	600
350-245	Erbe VIO international *)	ErgoLAP bipolar	600
351-040	BOWA, Erbe international *), Martin International, Valleylab, Conmed	BOWA forceps, EU flat connector	550
351-045	BOWA, Erbe international *), Martin International, Valleylab, Conmed	US forceps	550
351-051	BOWA, Erbe international *), Martin International, Valleylab, Conmed	BOWA forceps, EU flat connector	550
351-245	BOWA, Erbe international *), Martin International, Valleylab, Conmed	ErgoLAP bipolar	600
351-345	BOWA	ErgoLAP bipolar	600
351-445	BOWA, Erbe international *), Martin International, Valleylab, Conmed	Aesclapal forceps, flat connector	550
351-545	BOWA, Erbe international *), Martin International, Valleylab, Conmed	Aesclapal forceps, flat connector	550
352-145	BOWA	Resectoscope, Storz	1000
352-245	BOWA	Bipolar cable, Storz Laparoscopy, COMFORT	550
353-040	Erbe VIO international *)	BOWA forceps, EU flat connector	550
353-045	Erbe VIO international *)	US forceps	550
353-050	Erbe VIO international *)	EU forceps, EU flat connector	550
354-145	BOWA	Resectoscope, Wolf	1000
355-030	Valleylab	BOWA forceps, EU flat connector	550
355-031	Valleylab	BOWA forceps, EU flat connector	550
355-145	BOWA	Resectoscope, Olympus	1000
358-245	BOWA	ERGO 315R	200
378-045	BOWA, Erbe international *), Martin International, Valleylab, Conmed	Bipolar scissors (BIZZER, Ethicon)	300
378-145	Erbe VIO / ICC / ACC / T series	Bipolar scissors (BIZZER, Ethicon)	300
378-245	Martin	Bipolar scissors (BIZZER, Ethicon)	300
378-345	Martin	Storz bipolar laparoscopic handle	550
378-445	Martin	Aesclapal forceps, flat connector	550
Adapters			
340-000	Erbe	Bipolar cable (Valleylab/ Conmed) 2-pin	1000
901-210	BOWA ARC 250 / 303 / 350 / 400	For connecting up Erbe accessories	1000
901-220	BOWA ARC 250 / 303 / 350 / 400	For connecting up Martin accessories	1000

*) depending on the respective socket configuration

Contraindications:

Do not use HF cables / HF adapters in direct contact with the heart, central circulatory system or the central nervous system.

Do not use HF cables / HF adapters where its surgical techniques are contraindicated.

Do not use bipolar application if, in the opinion of an experienced physician or according to current professional literature, such use would cause endangerment of the patient, for example due to the general condition of the patient, or if other contraindications are present.

HF unit/power settings:

The maximum permissible voltage can be found in the table above.

2 Safety instructions

2.1 Device-related

- Goods which are delivered non-sterile must be cleaned and sterilised before use! Sterilisable up to 137°C.
- Do not sterilise the product in hot air.
- Check the device for damage after reprocessing and before use
- **CAUTION:** Visual inspection on its own is not sufficient for ensuring insulation integrity.
- If necessary, check the functional integrity with a suitable test device (e.g. 050-230).
- Do not use the product if it is damaged.
- Do not carry out any repair work or maintenance on the device.
- Ensure that the product is cleaned, disinfected and sterilised correctly.
- In order to avoid the risk of injuries and electric shocks for the patient and operating personnel, ensure that the power supply is switched off before connecting or disconnecting the adapter / HF cables and accessories to/from the electrosurgical device.
- Follow the instructions for use for the HF generator and related instruments.
- To disconnect, hold the HF cable by the plug and never pull directly on the HF cable itself.
- Risk of injury from sharp edges.

2.2 Use-related

- Users must be trained in, and familiar with, the basic principles & rules of the application, and risks of HF surgery.
- The products may only be used at medical facilities by trained medical personnel.
- Make sure that the cables do not come into contact with the patient.
- Position antistatic cloths on areas of the patient where substantial sweating is expected or there is skin-to-skin contact.
- Ensure that there are no explosive gases or inflammable fluids and materials in the surgical area.
- In the case of patients with pacemakers, consult a cardiologist before the application.
- Route the cable without spirals (formation of loops); keep it clear of other cables and not in contact with the patient.
- Comply with the recommended power settings and maximum voltage.
- Check that the cable is correctly connected to the instrument and the HF generator with the appropriate connector.
- Do not overload the HF cables beyond their mechanical limits (e.g., excessive bending or crushing as a result of being run over by an instrument trolley).
- Caution: When the AUTOSTART function is in use, HF application is activated by skin contact.
- Ensure that the cable length for bipolar applications in AUTOSTART mode does not exceed 4.5 m.
- To avoid bodily injuries resulting from inadvertent HF activation, do not place active HF instruments on the patient.
- The "BF"/"CF" application device of the HF device being used is extended by the connected instrument.
- Report any incidents or near misses to your country's medical devices regulatory agency. When doing so, follow your own facility's internal reporting system. Notify your local BOWA specialist retailer in such cases.

2.3 EMC information

Medical electrical devices are subject to special EMC precautions, so please follow the below-mentioned instructions.

- BOWA accessories are intended only for connection to the specified HF devices.
- Using the accessories with medical devices from other manufacturers can result in higher emission levels or reduced interference immunity.

i NOTE
Combining medical devices is only safe if <ul style="list-style-type: none"> ➤ the relevant instructions for use permit the desired combination or ➤ the intended use and the interface specification of the products used in the combination are compatible. Strict attention must be given to the instructions for use and interface specifications of medical devices used in combination.

3 Reprocessing

3.1 General information

The following information on the allowed number of reprocessing cycles should be regarded as a guideline. The actual number may vary depending on the demands placed on the product. Manual reprocessing methods are not recommended due to their significantly lower effectiveness. Follow the chemical manufacturers' statements on concentrations and working times. If other types of cleaning agents and disinfectants are used, the medical device manufacturer cannot be held liable.

3.2 Product-specific information

Do not clean cables and adapters in an ultrasonic bath.

Isolated parts must not come into contact with hard, sharp or heavy products, as this could compromise the electrical isolation.

Do not use wire brushes, abrasives or other products which could damage the surface.

i NOTE

The BOWA COMFORT instruments with instrument identification Plug'n Cut COMFORT are compatible with the following software versions of your generator:

- ARC 400 starting from V 2.2.1
- ARC 350 starting from V 2.1.1

If you operate your generator with an older software version, a software update will be required in order to use the BOWA COMFORT instruments with Plug'n Cut COMFORT. The use with software version 1.x.x is not possible because an update to a higher version is not supported.

European Medical Device Regulation
This product complies with Regulation 2017/745 on medical devices. **CE**

1 Intended Purpose

Electrosurgical equipment for cutting and coagulation of tissue.

Indications:

Bipolar HF-cables / HF-adapters form the electrical connection between a HF surgical unit and a bipolar application part to transmit HF current.

Reprocessing cycles	300 (Connection cables); 100 (COMFORT cables); 75 (Adapter); 20 (REF 352-145); 50 (REF 354-145, 355-145)
Reprocessing step	Description
Manual pre-cleaning	
1 Soaking¹	Soak the product at room temperature (<25°C) for at least 15 minutes immediately, but no later than 2 hours after use. Only use enzymatic, aldehyde-free disinfectants which are suitable for disinfecting medical devices (e.g., DGHM or FDA approval or CE marking). Remove all visible impurities with a soft plastic brush. BOWA recommends the use of Gigasept® Instru AF (Schülke & Mayr GmbH).
2 Rinsing	Rinse the product thoroughly for at least 1 minute at room temperature (<25°C) under running tap water (at least drinking water quality). Allow residual water to drip off sufficiently.
Automatic cleaning, thermal disinfection and drying	
3 Automatic cleaning²	Use a washer-disinfector (WD) with tested effectiveness (according to ISO 15883) and a neutral to alkaline (max pH 11.5) enzymatic cleaning agent. Depending on the concentration, agents including alcohol- and/or aldehyde-containing ingredients may be used. Ensure that the cables are not kinked or jammed. BOWA recommends the use of neodisher® MediClean forte (Chemische Fabrik Dr. Weigert) at 55°C for 10 minutes. The intermediate rinsing must include at least two rinsing steps with fully demineralised water (deionised water) (at least 1 minute each at >10°C).
4 Thermal disinfection³	An A ₀ value of >3000 must be maintained. BOWA recommends a temperature of 90°C for at least 5 minutes and the use of demineralised water.
5 Drying	Drying is carried out according to the WD programme and depends on the total demand. The maximum temperature is 100°C for 25 minutes. If necessary, dry with filtered compressed air at max. 3 bar.
Inspection	After machine cleaning and disinfection, carry out a visual inspection for residues. If necessary, repeat the reprocessing process.
Steam sterilisation	
6 Steam sterilisation⁴	BOWA recommends steam sterilisation using the fractionated vacuum method with a sterilisation time of 3–20 minutes and a temperature of 134–137°C in a suitable sterilisation packaging.

BOWA has validated:

1 Soaking with GigaseptR Instru AF (Schülke & Mayr GmbH), 3% v/v

2 Neodisher® MediClean forte (Chemische Fabrik Dr. Weigert), 0,5% v/v

3 Thermal disinfection at 90°C for 5 minutes

4 Half-cycle with 3 pre-vacuum phases at 132°C for 1.5 minutes, packaging: double packaging, in paper/film sterilisation bag

4 Environmental conditions

4.1 Transport and storage

BOWA recommends storing the products in their original packaging until their first use. Store the products in a location where they are protected against direct sunlight. Transport or store the instrument in a dry environment.

4.2 Operation

Use the product at room temperature, protected from mechanical force such as pulling or bending.

5 Disposal

Dispose of the medical devices, packaging material and accessories in accordance with the respectively applicable country-specific regulations and laws.

6 Symbols to bear in mind

Symbol	Meaning	Symbol	Meaning
	Follow the instructions for use		Keep away from direct sunlight
REF	Product number		Keep dry
LOT	Batch		Conforms to Russian standards
QTY	Quantity		Manufacturer
CE	CE mark	CE 0123	CE mark and identification number of Notified Body (class IIb)
R_{Only}	Caution: Sale and prescription restricted to physicians. For the treating physician only.		WEEE symbol: Indicates separate collection for waste of electric and electronic equipment.
MD	Medical device		Caution

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Copyright by BOWA-electronic GmbH & Co. KG, Gomaringen | Germany